

CLONING SMALL HAIRPINS INTO LENTIVIRAL VECTORS

Annealing Step: in a single tube, mix the following:

23 μ l ddH₂O
1 μ l sense oligo (100 μ M)
1 μ l antisense oligo (100 μ M)
25 μ l 2X Annealing Buffer (recipe see right)
50 μ l total reaction volume

2X Annealing Buffer: 200 mM potassium acetate 60 mM HEPES-KOH pH 7.4 4 mM Mg-acetate
--

Incubate the reaction for 5 minutes at 95°C and let the oligos gradually anneal to RT. These can then be placed on ice (or stored at -20°C) until used in ligation step. *This is the annealed oligo stock; prior to use, it will be necessary to prepare a 1:20 dilution.*

Ligation Step: in a single tube, mix in the following order

x μ l ddH₂O
50-100 ng digested, gel purified vector
1 μ l of 1:20 diluted annealed oligos
1 μ l 10X T4 DNA ligase buffer
1 μ l T4 DNA ligase
10 μ l total reaction volume

Incubate the above reaction at room temperature for 1-4 hours and transform 2.5 μ l of the ligation reaction into DH5 α cells.